



Innovative Infusion Systems, Since 1994



Applicable
syringe:
5, 10, 20, 30,
50/60 ml

0.01 ml/h
increments for
flow rate,
bolus rate,
KVO rate



Syringe Pump

KL-605T



PROMEDIC.GE
DISTRIBUIDOR OFICIAL

Beijing KellyMed Co., Ltd.

Syringe Pump

KL-605T



Specifications

Model	KL-605T
Syringe Size	5, 10, 20, 30, 50/60 ml
Applicable Syringe	Compatible with syringe of any standard
VTBI	1-1000 ml (in 0.1, 1, 10 ml increments)
Flow Rate	Syringe 5 ml: 0.1-100 ml/h (in 0.01, 0.1, 1, 10 ml/h increments) Syringe 10 ml: 0.1-300 ml/h Syringe 20 ml: 0.1-600 ml/h Syringe 30 ml: 0.1-800 ml/h Syringe 50/60 ml: 0.1-1200 ml/h
Bolus Rate	5 ml: 0.1-100 ml/h (in 0.01, 0.1, 1, 10 ml/h increments) 10 ml: 0.1-300 ml/h 20 ml: 0.1-600 ml/h 30 ml: 0.1-800 ml/h 50/60 ml: 0.1-1200 ml/h
Anti-Bolus	Automatic
Accuracy	± 2% (mechanical accuracy ≤ 1%)
Infusion Mode	Flow rate: ml/min, ml/h Time-based Body weight
KVO Rate	0.1-1 ml/h (in 0.01 ml/h increments)
Alarms	Occlusion, near empty, end program, low battery, end battery, AC power off, motor malfunction, system malfunction, reminder alarm, pressure sensor error, syringe installation error, syringe drop off
Additional Features	Real-time infused volume, automatic power switching, automatic syringe identification, mute button, purge, bolus, anti-bolus, system memory, history log
Drug Library	Available
Occlusion Sensitivity	High, medium, low
Power Supply, AC	110/230 V, 50/60 Hz, 20 VA
Battery	14.8 V, rechargeable
Battery Life	8 hours at 5 ml/h
Working Temperature	5-40°C
Relative Humidity	20-90%
Atmospheric Pressure	700-1060 hpa
Size	245*120*115 mm
Weight	2.5 kg
Safety Classification	Class II, type BF

Beijing KellyMed Co., Ltd.

Add: 6R International Metro Center, No. 3 Shilipu,
Chaoyang District, Beijing, 100025, China

Tel: +86-10-82490385

Fax: +86-10-65587908

E-mail: export@kelly-med.com

Web: www.kelly-med.com

Distributor: